Presentation on. 26-04-2025

EfficientNet-Based and YOLO-Driven **Brain Tumor Detection and** Segmentation

Authors:

Najla Musthafa¹, Mohammed Aflah², Minhaj Akavalappil², Mohammed Jasim A², Mohammed Aseel², Shanid Malayil³, Mubeena A K¹.

¹Assistant Professor, Dept. of CSE, MEA Engineering College, Kerala, India

²Dept. of CSE, MEA Engineering College, Kerala, India

³Associate Professor, Dept. of CSE, MEA Engineering College, Kerala, India















- Introduction
- Background Information
- Proposed System
- Methodology
- Results
- Discussion
- Conclusion
- References

Introduction



- Brain tumors are life-threatening and require early, accurate diagnosis to improve survival rates.
- Manual MRI diagnosis is time-consuming, subjective, and prone to human error.
- Al and Deep Learning offer automation, speed, and increased diagnostic accuracy.
- EfficientNet-B0 is used to classify tumors into four categories (glioma, meningioma, pituitary, no tumor).
- YOLOv8 performs real-time segmentation to localize the tumor in MRI scans.
- This integrated system supports radiologists with fast and reliable second opinions.




Background Information

- MRI is the standard imaging tool for brain tumors due to its detailed visualization of soft tissues.
- Traditional ML methods (like SVMs, kNN) rely on handcrafted features and struggle with complex medical data.
- CNNs automatically extract deep features and improve image classification and segmentation performance.
- Transfer Learning enables pre-trained models to work effectively on small medical datasets.
- EfficientNet balances model size and accuracy through compound scaling.
- YOLOv8 ensures fast, accurate segmentation—ideal for clinical use and real-time diagnosis.

Proposed System







Methodology

1. Dataset Collection:

Collect **5000 MRI images** across 4 categories — glioma tumor, meningioma tumor, pituitary tumor, and no tumor.

2. Data Preprocessing:

- **Augmentation:** Apply rotation, flipping, zooming for data diversity.
- Annotation: Tumor regions manually marked for segmentation.

3. Model Selection:

- Use **EfficientNet-B0** (pretrained on ImageNet) for classification.
- Use YOLOv8 for segmentation (to detect and outline tumor regions).

4. Training:

- EfficientNet-B0: Fine-tuned using Transfer Learning for 4-class classification.
- YOLOv8: Trained to detect and segment tumors using annotated images.





5. Loss Functions:

Classification Loss: Categorical Cross-Entropy:

$$L = -\sum (y \text{ true } * \log(y \text{ pred}))$$

Segmentation Loss: Intersection over Union (IoU):

$$IoU = |A \cap B| / |A \cup B|$$

6. Evaluation:

- Classification Metrics: Accuracy, Sensitivity, Specificity, F1-score.
- **Segmentation Metric:** IoU (measures overlap between predicted and actual tumor area).

7. Deployment:

 Plan to develop a simple application for radiologists to upload MRI scans and get instant tumor detection and classification results.



Methodology - Classification

- Images resized to 150×150 for EfficientNet input.
- Transfer learning applied using weights from ImageNet.
- Final layer: **Softmax** for 4-class probability distribution.

$$\hat{y}_i = rac{e^{z_i}}{\sum_j e^{z_j}}$$

- Each class probability = (exponentiated value) / (sum of all exponentiated values).
- \circ z_i = raw score for class i.
- **Example:** If scores for 4 classes are: [2, 1.5, 0.5, -0.5]

$$\hat{y}_{glioma} = rac{e^2}{e^2 + e^{1.5} + e^{0.5} + e^{-0.5}}$$
 This

This gives probabilities like [0.7, 0.2, 0.08, 0.02].

Loss Function: Categorical Cross-Entropy

$$L = -\sum(y_true * log(y_pred))$$



Methodology - Segmentation


- Images resized to 416×416 for YOLOv8.
- Annotated tumor regions (bounding boxes and masks) fed into model.
- Segmentation performance measured using: **IoU** (Intersection over Union):

$$IoU = |A \cap B| / |A \cup B|$$

- A: The predicted tumor region.
- B: The actual tumor region (from the dataset).
- YOLOv8 provides real-time segmentation outputs.



Confusion Matrix



Result

Accurate Tumor Classification

- EfficientNet-B0 achieves final training accuracy of approximately 99.7%.
- Validation accuracy reaches 98.3%, indicating excellent generalization to unseen data.
- High **F1-score of 96%** across all classes confirms **balanced performance**.


Efficient Segmentation

- YOLOv8 segmentation model achieves a mean Average Precision (mAP) of 81.2%.
- Tumor regions (glioma, meningioma, pituitary) are accurately highlighted in MRI scans.
- Real-time segmentation ensures faster processing.



Result

- Best Model: EfficientNetB0 with Adam.
- Why EfficientNetB0?
 - Higher accuracy than MobileNetV2, VGG16, and ResNet50.
 - More efficient with fewer parameters than ResNet50 and VGG16.
 - Better feature extraction with compound scaling.
- Key Takeaway: EfficientNetB0-Adam is the best choice for brain tumor detection, ensuring high accuracy with efficiency.





Result

No Tumor



Classification



The Model predicts that it is a Glioma Tumor

Segmentation







Discussion And Future Scope

Diagnosis & Surgery

Automates tumor detection and aids in surgical planning with precise segmentation.

Remote Access

Enables diagnosis via cloud platforms and telemedicine in underserved areas.

Mobile Health

Integrates into smartphone apps for quick, accessible tumor screening.

Emergency Use

Real-time analysis supports fast decisions in critical care settings.

Monitoring & Research

Tracks treatment progress and supports medical AI research.

Explainability & Expansion

Adds interpretability for clinical trust and combines data from CT, PET, and genetics for better accuracy.



Conclusion

- The proposed system successfully combines EfficientNet-B0 for classification and YOLOv8 for segmentation of brain tumors.
- It achieves high classification accuracy (99.7%) and strong segmentation performance (81.2% mAP).
- Transfer Learning significantly reduces training time and enhances performance on small datasets.
- The framework delivers **real-time results**, aiding radiologists in faster and more accurate diagnoses.
- The system supports automated medical diagnostics, reducing radiologist workload and minimizing human error.
- This research highlights the potential of Al-driven healthcare solutions to improve early detection and patient outcomes.

